





Features of the Authentic EES Product PROXIMATE® PPH Procedure for Prolapse and Hemorrhoids set	Features of the Counterfeit Product
A. Authentic EES Product has a bar code on one side of the flaps. Note: Both flaps must be inspected. (See Picture 1)	A. Counterfeit Product does not have a bar code on the flap.
 <p>Picture 1</p>	
Authentic EES Product: One side of the flaps has a bar code	Counterfeit Product: No bar code on both flaps
B. Authentic EES Product has a batch number on the firing trigger. The batch number should be visible in an unopened blister by removing the blister from the carton. (see picture 2)	B. Counterfeit Product does not have a batch number on the firing trigger.
 <p>Picture 2</p>	
Authentic EES Product: Batch number on firing trigger	Counterfeit Product: No batch number on firing trigger

